

# Applications



## Gingiva UV 2.0

- ✓ Elastic and tear-resistant
- ✓ Dimensional stability during print
- ✓ Limited shrinkage
- ✓ Aging resistance
- ✓ Long-term storage
- ✓ Odorless
- ✓ Natural tone



Gingiva UV 2.0



Gingiva UV 2.0 has excellent dimensional stability, elasticity, and tear-resistance for 3D reproduction of functional gingival model segments in a digital workflow in conjunction with our Model UV materials. It remains elastic after extended use without shrinkage or degradation and is odorless after the printing process.

**Color** Pink ●

**Specification** 500g Bottle / 1000g Bottle

**Certifications** NMPA Free / CE Free / FDA

**Machine Compatibility**



A2D



A2D 4K

Properties	Units	Standard	Result
Hardness	Shore A	ASTM D2240	38
Elongation at Break	%	ASTM D638	134.6–143.1
Tear Resistance	N/mm	ASTM D624	7.8–10.5